

## Main specification

4,3" LCD display	yes
USB port for easy upload of software/download of logged data	yes
Extensive integrated logging capabilities with access to detailed log directly on the display	yes
Monitoring of multiple alarm conditions	yes
Integrated memory for 30 years	yes
Available in multiple languages	yes
Battery backup for alarms, logging and temperature display in case of power cut	yes
Auto cycle in case of probe failure	yes
3 level of password protection	yes

## Optional features

Prepared for connection of 2 additional probes	yes
Prepared for electronic lock with access via RFID card	yes
Prepared RS485 via Modbus	yes
Prepared for connection of GSM alarm module	yes

## Alarm features

Visual and acoustic alarm	yes
Adjustable high/low temperature alarm	yes
Power failure alarm	yes
Probe failure alarm	yes
Instrument failure alarm	yes
Open door alarm	yes
Contact for remote alarm	yes
Shows all alarms	yes