

Box size 32

Model	Box size (height mm)	Rack model & total quantity	No. of boxes per rack (depth x height)	Total no. of boxes	Exterior dimensions of rack (wxdxh mm)
ULUF P10MV	32	-	-	9	-
ULUF P50MV	32	2 x USR34-38	3 x 8	48	139 x 414 x 292
ULUF P90MV	32	4 x USR34-38	3 x 8	96	139 x 414 x 292
ULUF P400MV	32	6 x UDR34-47	4 x 7	168	140 x 567 x 270
		6 x UDR34-49	4 x 9	216	140 x 567 x 347
		Total: 12 racks	-	(=384)	-
ULUF P510MV	32	8 x UDR34-47	4 x 7	224	140 x 567 x 270
		8 x UDR34-49	4 x 9	288	140 x 567 x 347
		Total: 16 racks	-	(=512)	-
ULUF P620MV	32	10 x UDR34-47	4 x 7	280	140 x 567 x 270
		10 x UDR34-49	4 x 9	360	140 x 567 x 347
		Total: 20 racks	-	(=640)	-
ULUF P830MV	32	14 x UDR34-47	4 x 7	392	140 x 567 x 270
		14 x UDR34-49	4 x 9	504	140 x 567 x 347
		Total: 28 racks	-	(=896)	-
ULUF 450-2M	32	15 x UDR34-46	4 x 6	360	140 x 567 x 232
ULUF 750-2M	32	25 x UDR34-46	4 x 6	600	140 x 567 x 232

Box size 40

Model	Box size (height mm)	Rack model & total quantity	No. of boxes per rack (depth x height)	Total no. of boxes	Exterior dimensions of rack (wxdxh mm)
ULUF P10MV	40	-	-	7	-
ULUF P50MV	40	2 x UDR40-36	3 x 6	36	140 x 434 x 268
ULUF P90MV	40	4 x UDR40-36	3 x 6	72	140 x 434 x 268
ULUF P400MV	40	6 x UDR40-46	4 x 6	144	140 x 567 x 268
		6 x UDR40-47	4 x 7	168	140 x 567 x 312
		Total: 12 racks	-	(=312)	-
ULUF P510MV	40	8 x UDR40-46	4 x 6	192	140 x 567 x 268
		8 x UDR40-47	4 x 7	224	140 x 567 x 312
		Total: 16 racks	-	(=416)	-
ULUF P620MV	40	10 x UDR40-46	4 x 6	240	140 x 567 x 268
		10 x UDR40-48	4 x 8	320	140 x 567 x 357
		Total: 20 racks	-	(=560)	-
ULUF P830MV	40	14 x UDR40-46	4 x 6	336	140 x 567 x 268
		14 x UDR40-48	4 x 8	448	140 x 567 x 357
		Total: 28 racks	-	(=784)	-
ULUF 450-2M	40	15 x UDR40-45	4 x 5	300	140 x 567 x 223
ULUF 750-2M	40	25 x UDR40-45	4 x 5	500	140 x 567 x 223

Box size 50

Model	Box size (height mm)	Rack model & total quantity	No. of boxes per rack (depth x height)	Total no. of boxes	Exterior dimensions of rack (wxdxh mm)
ULUF P10MV	50	-	-	5	-
ULUF P50MV	50	2 x UDR50-35	3 x 5	30	140 x 432 x 280
ULUF P90MV	50	4 x UDR50-35	3 x 5	60	140 x 432 x 280
ULUF P400MV	50	6 x UDR50-45	4 x 5	120	140 x 570 x 280
		6 x UDR50-46	4 x 6	144	140 x 570 x 336
		Total: 12 racks	-	(=264)	-
ULUF P510MV	50	8 x UDR50-45	4 x 5	160	140 x 570 x 280
		8 x UDR50-46	4 x 6	192	140 x 570 x 336
		Total: 16 racks	-	(=352)	-
ULUF P620MV	50	10 x UDR50-45	4 x 5	200	140 x 570 x 280
		10 x UDR50-46	4 x 6	240	140 x 570 x 336
		Total: 20 racks	-	(=440)	-
ULUF P830MV	50	14 x UDR50-45	4 x 5	280	140 x 570 x 280
		14 x UDR50-46	4 x 6	336	140 x 570 x 336
		Total: 28 racks	-	(=616)	-
ULUF 450-2M	50	15 x UDR50-44	4 x 4	240	140 x 570 x 224
ULUF 750-2M	50	25 x UDR50-44	4 x 4	400	140 x 570 x 224

Box size 75

Model	Box size (height mm)	Rack model & total quantity	No. of boxes per rack (depth x height)	Total no. of boxes	Exterior dimensions of rack (wx dxh mm)
ULUF P10MV	75	-	-	4	-
ULUF P50MV	75	2 x UDR75-33	3 x 3	18	140 x 420 x 242
ULUF P90MV	75	4 x UDR75-33	3 x 3	36	140 x 420 x 242
ULUF P400MV	75	6 x UDR75-43	4 x 3	72	140 x 558 x 240
		6 x UDR75-44	4 x 4	96	140 x 558 x 335
		Total: 12 racks		(=168)	-
ULUF P510MV	75	8 x UDR75-43	4 x 3	96	140 x 558 x 240
		8 x UDR75-44	4 x 4	128	140 x 558 x 335
		Total: 16 racks	-	(=224)	-
ULUF P620MV	75	10 x UDR75-43	4 x 3	120	140 x 558 x 240
		10 x UDR75-44	4 x 4	160	140 x 558 x 335
		Total: 20 racks	-	(=280)	-
ULUF P830MV	75	14 x UDR75-43	4 x 3	168	140 x 558 x 240
		14 x UDR75-44	4 x 4	224	140 x 558 x 335
		Total: 28 racks	-	(=392)	-
ULUF 450-2M	75	15 x UDR75-43	4 x 3	180	140 x 558 x 240
ULUF 750-2M	75	25 x UDR75-42	4 x 2	200	140 x 558 x 162

Box size 80

Model	Box size (height mm)	Rack model & total quantity	No. of boxes per rack (depth x height)	Total no. of boxes	Exterior dimensions of rack (wx dx h mm)
ULUF P10MV	80	-	-	3	-
ULUF P50MV	80	2 x UDR80-23	2 x 3	12	140 x 301 x 254
ULUF P90MV	80	4 x UDR80-23	2 x 3	24	140 x 301 x 254
ULUF P400MV	80	6 x UDR80-43	4 x 3	72	140 x 567 x 254
		6 x UDR80-44	4 x 4	96	140 x 567 x 339
		Total: 12 racks	-	(=168)	-
ULUF P510MV	80	8 x UDR80-43	4 x 3	96	140 x 567 x 254
		8 x UDR80-44	4 x 4	128	140 x 567 x 339
		Total: 16 racks	-	(=224)	-
ULUF P620MV	80	10 x UDR80-43	4 x 3	120	140 x 567 x 254
		10 x UDR80-44	4 x 4	160	140 x 567 x 339
		Total: 20 racks	-	(=280)	-
ULUF P830MV	80	14 x UDR80-43	4 x 3	168	140 x 567 x 254
		14 x UDR80-44	4 x 4	224	140 x 567 x 339
		Total: 28 racks	-	(=392)	-
ULUF 450-2M	80	15 x UDR80-43	4 x 3	180	140 x 567 x 254
ULUF 750-2M	80	25 x UDR80-43	4 x 3	200	140 x 567 x 254

Box size 90/96/100

Model	Box size (height mm)	Rack model & total quantity	No. of boxes per rack (depth x height)	Total no. of boxes	Exterior dimensions of rack (wxdxh mm)
ULUF P10MV	90/96/100	-	-	3	-
ULUF P50MV	90/96/100	2 x UDR100-32	3 x 2	12	140 x 420 x 212
ULUF P90MV	90/96/100	4 x UDR100-32	3 x 2	24	140 x 420 x 212
ULUF P400MV	90/96/100	6 x UDR100-42	4 x 2	48	140 x 558 x 212
		6 x UDR100-43	4 x 3	72	140 x 558 x 318
		Total: 12 racks	-	(=120)	-
ULUF P510MV	90/96/100	8 x UDR100-42	4 x 2	64	140 x 558 x 212
		8 x UDR100-43	4 x 3	96	140 x 558 x 318
		Total: 16 racks	-	(=160)	-
ULUF P620MV	90/96/100	10 x UDR100-42	4 x 2	80	140 x 558 x 212
		10 x UDR100-43	4 x 3	120	140 x 558 x 318
		Total: 20 racks	-	(=200)	-
ULUF P830MV	90/96/100	14 x UDR100-42	4 x 2	112	140 x 558 x 212
		14 x UDR100-43	4 x 3	168	140 x 558 x 318
		Total: 28 racks	-	(=280)	-
ULUF 450-2M	90/96/100	15 x UDR100-42	4 x 2	120	140 x 558 x 212
ULUF 750-2M	90/96/100	25 x UDR100-42	4 x 2	200	140 x 558 x 212