

Box size 32

Model	Box size (height mm)	Rack model & total quantity	No. of boxes per rack (depth x height)	Total no. of boxes	Exterior dimensions of rack (wxdxh mm)
<b>MUF 40</b>	32	2 x UDR34-27	2 x 7	28	140 x 301 x 270
<b>ULUF P10</b>	32	-	-	9	-
<b>ULUF P50</b>	32	2 x USR34-38	3 x 8	48	139 x 414 x 292
<b>ULUF P90</b>	32	4 x USR34-38	3 x 8	96	139 x 414 x 292
<b>ULUF P390</b>	32	9 x UDR34-47	4 x 7	252	140 x 567 x 270
		3 x UDR34-49	4 x 9	108	140 x 567 x 347
		Total: 12 racks	-	(=360)	-
<b>ULUF P390 GG</b>	32	6 x UDR34-47	4 x 7	168	140 x 567 x 270
		6 x UDR34-49	4 x 9	216	140 x 567 x 347
		Total: 12 racks	-	(=384)	-
<b>ULUF P500</b>	32	12 x UDR34-47	4 x 7	336	140 x 567 x 270
		4 x UDR34-49	4 x 9	144	140 x 567 x 347
		Total: 16 racks	-	(=480)	-
<b>ULUF P500 GG</b>	32	8 x UDR34-47	4 x 7	224	140 x 567 x 270
		8 x UDR34-49	4 x 9	288	140 x 567 x 347
		Total: 16 racks	-	(=512)	-
<b>ULUF P610</b>	32	15 x UDR34-47	4 x 7	420	140 x 567 x 270
		5 x UDR34-49	4 x 9	180	140 x 567 x 347
		Total: 20 racks	-	(=600)	-
<b>ULUF P610 GG</b>	32	10 x UDR34-47	4 x 7	280	140 x 567 x 270
		10 x UDR34-49	4 x 9	360	140 x 567 x 347
		Total: 20 racks	-	(=640)	-
<b>ULUF P820</b>	32	21 x UDR34-47	4 x 7	588	140 x 567 x 270
		7 x UDR34-49	4 x 9	252	140 x 567 x 347
		Total: 28 racks	-	(=840)	-
<b>ULUF P820 GG</b>	32	14 x UDR34-47	4 x 7	392	140 x 567 x 270
		14 x UDR34-49	4 x 9	504	140 x 567 x 347
		Total: 28 racks	-	(=896)	-
<b>ULUF 15</b>	32	-	-	9	-
<b>ULUF 65</b>	32	2 x UDR34-38	3 x 8	48	140 x 434 x 309
<b>ULUF 60</b>	32	2 x UDR34-38	3 x 8	48	140 x 434 x 309
<b>ULUF 120</b>	32	4 x UDR34-36	3 x 6	72	140 x 434 x 232
<b>ULUF 125</b>	32	4 x UDR34-36	3 x 6	72	140 x 434 x 232
<b>ULUF 400</b>	32	12 x UDR34-46	4 x 6	288	140 x 567 x 232
		3 x UDR34-48	4 x 8	96	140 x 567 x 309
		Total: 15 racks	-	(=384)	-
<b>ULUF 450-2M</b>	32	15 x UDR34-46	4 x 6	360	140 x 567 x 232
<b>ULUF 500</b>	32	16 x UDR34-46	4 x 6	384	140 x 567 x 232
		4 x UDR34-48	4 x 8	128	140 x 567 x 309
		Total: 20 racks	-	(=512)	-
<b>ULUF 550-2M</b>	32	20 x UDR34-46	4 x 6	480	140 x 567 x 232
<b>ULUF 700</b>	32	20 x UDR34-46	4 x 6	480	140 x 567 x 232
		5 x UDR34-48	4 x 8	160	140 x 567 x 309
		Total: 25 racks	-	(= 640)	-
<b>ULUF 750-2M</b>	32	25 x UDR34-46	4 x 6	600	140 x 567 x 232

Box size 40

Model	Box size (height mm)	Rack model & total quantity	No. of boxes per rack (depth x height)	Total no. of boxes	Exterior dimensions of rack (wxdxh mm)
<b>MUF 40</b>	40	2 x UDR40-26	2 x 7	28	140 x 301 x 268
<b>ULUF P10</b>	40	-	-	7	-
<b>ULUF P50</b>	40	2 x UDR40-36	3 x 6	36	140 x 434 x 268
<b>ULUF P90</b>	40	4 x UDR40-36	3 x 6	72	140 x 434 x 268
<b>ULUF P390</b>	40	9 x UDR40-46	4 x 6	216	140 x 567 x 268
		3 x UDR40-48	4 x 8	96	140 x 567 x 357
		Total: 12 racks	-	(=312)	-
<b>ULUF P390 GG</b>	40	6 x UDR40-46	4 x 6	144	140 x 567 x 268
		6 x UDR40-47	4 x 7	168	140 x 567 x 312
		Total: 12 racks	-	(=312)	-
<b>ULUF P500</b>	40	12 x UDR40-46	4 x 6	288	140 x 567 x 268
		4 x UDR40-48	4 x 8	128	140 x 567 x 357
		Total: 16 racks	-	(=416)	-
<b>ULUF P500 GG</b>	40	8 x UDR40-46	4 x 6	192	140 x 567 x 268
		8 x UDR40-47	4 x 7	224	140 x 567 x 312
		Total: 16 racks	-	(=416)	-
<b>ULUF P610</b>	40	15 x UDR40-46	4 x 6	360	140 x 567 x 268
		5 x UDR40-48	4 x 8	160	140 x 567 x 268
		Total: 20 racks	-	(=520)	-
<b>ULUF P610 GG</b>	40	10 x UDR40-46	4 x 6	240	140 x 567 x 268
		10 x UDR40-48	4 x 8	320	140 x 567 x 357
		Total: 20 racks	-	(=560)	-
<b>ULUF P820</b>	40	21 x UDR40-46	4 x 6	504	140 x 567 x 268
		7 x UDR40-48	4 x 8	224	140 x 567 x 357
		Total: 28 racks	-	(=728)	-
<b>ULUF P820 GG</b>	40	14 x UDR40-46	4 x 6	336	140 x 567 x 268
		14 x UDR40-48	4 x 8	448	140 x 567 x 357
		Total: 28 racks	-	(=784)	-
<b>ULUF 15</b>	40	-	-	7	-
<b>ULUF 65</b>	40	2 x UDR40-36	3 x 6	36	140 x 434 x 268
<b>ULUF 60</b>	40	2 x UDR40-36	3 x 6	36	140 x 434 x 268
<b>ULUF 120</b>	40	4 x UDR40-35	3 x 5	60	140 x 434 x 223
<b>ULUF 125</b>	40	4 x UDR40-35	3 x 5	60	140 x 434 x 223
<b>ULUF 400</b>	40	12 x UDR40-45	4 x 5	240	140 x 567 x 223
		3 x UDR40-47	4 x 7	84	140 x 567 x 312
		Total: 15 racks	-	(=324)	-
<b>ULUF 450-2M</b>	40	15 x UDR40-45	4 x 5	300	140 x 567 x 223
<b>ULUF 500</b>	40	16 x UDR40-45	4 x 5	320	140 x 567 x 223
		4 x UDR40-47	4 x 7	112	140 x 567 x 312
		Total: 20 racks	-	(=432)	-
<b>ULUF 550-2M</b>	40	20 x UDR40-45	4 x 5	400	140 x 567 x 223
<b>ULUF 700</b>	40	20 x UDR40-45	4 x 5	400	140 x 567 x 223
		5 x UDR40-47	4 x 7	140	140 x 567 x 312
		Total: 25 racks	-	(=540)	-
<b>ULUF 750-2M</b>	40	25 x UDR40-45	4 x 5	500	140 x 567 x 223

Box size 50

Model	Box size (height mm)	Rack model & total quantity	No. of boxes per rack (depth x height)	Total no. of boxes	Exterior dimensions of rack (wxdxh mm)
<b>MUF 40</b>	50	2 x UDR50-25	2 x 5	20	140 x 294 x 280
<b>ULUF P10</b>	50	-	-	5	-
<b>ULUF P50</b>	50	2 x UDR50-35	3 x 5	30	140 x 432 x 280
<b>ULUF P90</b>	50	4 x UDR50-35	3 x 5	60	140 x 432 x 280
<b>ULUF P390</b>	50	9 x UDR50-45	4 x 5	180	140 x 570 x 280
		3 x UDR50-46	4 x 6	72	140 x 570 x 336
		Total: 12 racks	-	(=252)	-
<b>ULUF P390 GG</b>	50	6 x UDR50-45	4 x 5	120	140 x 570 x 280
		6 x UDR50-46	4 x 6	144	140 x 570 x 336
		Total: 12 racks	-	(=264)	-
<b>ULUF P500</b>	50	12 x UDR50-45	4 x 5	240	140 x 570 x 280
		4 x UDR50-46	4 x 6	96	140 x 570 x 336
		Total: 16 racks	-	(=336)	-
<b>ULUF P500 GG</b>	50	8 x UDR50-45	4 x 5	160	140 x 570 x 280
		8 x UDR50-46	4 x 6	192	140 x 570 x 336
		Total: 16 racks	-	(=352)	-
<b>ULUF P610</b>	50	15 x UDR50-45	4 x 5	300	140 x 570 x 280
		5 x UDR50-46	4 x 6	129	140 x 570 x 336
		Total: 20 racks	-	(=420)	-
<b>ULUF P610 GG</b>	50	10 x UDR50-45	4 x 5	200	140 x 570 x 280
		10 x UDR50-46	4 x 6	240	140 x 570 x 336
		Total: 20 racks	-	(=440)	-
<b>ULUF P820</b>	50	21 x UDR50-45	4 x 5	420	140 x 570 x 280
		7 x UDR50-46	4 x 6	168	140 x 570 x 336
		Total: 28 racks	-	(=588)	-
<b>ULUF P820 GG</b>	50	14 x UDR50-45	4 x 5	280	140 x 570 x 280
		14 x UDR50-46	4 x 6	336	140 x 570 x 336
		Total: 28 racks	-	(=616)	-
<b>ULUF 15</b>	50	-	-	5	-
<b>ULUF 65</b>	50	2 x UDR50-35	3 x 5	30	140 x 432 x 280
<b>ULUF 60</b>	50	2 x UDR50-35	3 x 5	30	140 x 432 x 280
<b>ULUF 120</b>	50	4 x UDR50-34	3 x 4	48	140 x 420 x 228
<b>ULUF 125</b>	50	4 x UDR50-34	3 x 4	48	140 x 420 x 228
<b>ULUF 400</b>	50	12 x UDR50-44	4 x 4	192	140 x 570 x 224
		3 x UDR50-45	4 x 5	60	140 x 570 x 280
		Total: 15 racks	-	(=252)	-
<b>ULUF 450-2M</b>	50	15 x UDR50-44	4 x 4	240	140 x 570 x 224
<b>ULUF 500</b>	50	16 x UDR50-44	4 x 4	256	140 x 570 x 224
		4 x UDR50-45	4 x 5	80	140 x 570 x 280
		Total: 20 racks	-	(=336)	-
<b>ULUF 550-2M</b>	50	20 x UDR50-44	4 x 4	320	140 x 570 x 224
<b>ULUF 700</b>	50	20 x UDR50-44	4 x 4	320	140 x 570 x 224
		5 x UDR50-45	4 x 5	100	140 x 570 x 280
		Total: 25 racks	-	(=420)	-
<b>ULUF 750-2M</b>	50	25 x UDR50-44	4 x 4	400	140 x 570 x 224

Box size 75

Model	Box size (height mm)	Rack model & total quantity	No. of boxes per rack (depth x height)	Total no. of boxes	Exterior dimensions of rack (wxdxh mm)
<b>MUF 40</b>	75	2 x UDR75-23	2 x 3	12	140 x 294 x 240
<b>ULUF P10</b>	75	-	-	4	-
<b>ULUF P50</b>	75	2 x UDR75-33	3 x 3	18	140 x 420 x 242
<b>ULUF P90</b>	75	4 x UDR75-33	3 x 3	36	140 x 420 x 242
<b>ULUF P390</b>	75	9 x UDR75-43	4 x 3	108	140 x 558 x 240
		3 x UDR75-44	4 x 4	48	140 x 558 x 335
		Total: 12 racks	-	(=156)	-
<b>ULUF P390 GG</b>	75	6 x UDR75-43	4 x 3	72	140 x 558 x 240
		6 x UDR75-44	4 x 4	96	140 x 558 x 335
		Total: 12 racks	-	(=168)	-
<b>ULUF P500</b>	75	12 x UDR75-43	4 x 3	144	140 x 558 x 240
		4 x UDR75-44	4 x 4	64	140 x 558 x 335
		Total: 16 racks	-	(=208)	-
<b>ULUF P500 GG</b>	75	8 x UDR75-43	4 x 3	96	140 x 558 x 240
		8 x UDR75-44	4 x 4	128	140 x 558 x 335
		Total: 16 racks	-	(=224)	-
<b>ULUF P610</b>	75	15 x UDR75-43	4 x 3	180	140 x 558 x 240
		5 x UDR75-44	4 x 4	80	140 x 558 x 335
		Total: 20 racks	-	(=260)	-
<b>ULUF P610 GG</b>	75	10 x UDR75-43	4 x 3	120	140 x 558 x 240
		10 x UDR75-44	4 x 4	160	140 x 558 x 335
		Total: 20 racks	-	(=280)	-
<b>ULUF P820</b>	75	21 x UDR75-43	4 x 3	252	140 x 558 x 240
		7 x UDR75-44	4 x 4	112	140 x 558 x 335
		Total: 28 racks	-	(=364)	-
<b>ULUF P820 GG</b>	75	14 x UDR75-43	4 x 3	168	140 x 558 x 240
		14 x UDR75-44	4 x 4	224	140 x 558 x 335
		Total: 28 racks	-	(=392)	-
<b>ULUF 15</b>	75	-	-	4	-
<b>ULUF 65</b>	75	2 x UDR75-33	3 x 3	18	140 x 420 x 242
<b>ULUF 60</b>	75	2 x UDR75-33	3 x 3	18	140 x 420 x 242
<b>ULUF 120</b>	75	4 x UDR75-33	3 x 3	36	140 x 420 x 242
<b>ULUF 125</b>	75	4 x UDR75-33	3 x 3	36	140 x 420 x 242
<b>ULUF 400</b>	75	15 x UDR75-43	4 x 3	180	140 x 558 x 240
<b>ULUF 450-2M</b>	75	15 x UDR75-43	4 x 3	180	140 x 558 x 240
<b>ULUF 500</b>	75	20 x UDR75-43	4 x 3	240	140 x 558 x 240
<b>ULUF 550-2M</b>	75	20 x UDR75-43	4 x 3	240	140 x 558 x 240
<b>ULUF 700</b>	75	25 x UDR75-43	4 x 3	300	140 x 558 x 240
<b>ULUF 750-2M</b>	75	25 x UDR75-42	4 x 2	200	140 x 558 x 162

Box size 80

Model	Box size (height mm)	Rack model & total quantity	No. of boxes per rack (depth x height)	Total no. of boxes	Exterior dimensions of rack (wxdxh mm)
<b>MUF 40</b>	80	2 x UDR80-23	2 x 3	12	140 x 301 x 254
<b>ULUF P10</b>	80	-	-	3	-
<b>ULUF P50</b>	80	2 x UDR80-23	2 x 3	12	140 x 301 x 254
<b>ULUF P90</b>	80	4 x UDR80-23	2 x 3	24	140 x 301 x 254
<b>ULUF P390</b>	80	9 x UDR80-43	4 x 3	108	140 x 567 x 254
		3 x UDR80-44	4 x 4	48	140 x 567 x 339
		Total: 12 racks	-	(=156)	-
<b>ULUF P390 GG</b>	80	6 x UDR80-43	4 x 3	72	140 x 567 x 254
		6 x UDR80-44	4 x 4	96	140 x 567 x 339
		Total: 12 racks	-	(=168)	-
<b>ULUF P500</b>	80	12 x UDR80-43	4 x 3	144	140 x 567 x 254
		4 x UDR80-44	4 x 4	64	140 x 567 x 339
		Total: 16 racks	-	(=208)	-
<b>ULUF P500 GG</b>	80	8 x UDR80-43	4 x 3	96	140 x 567 x 254
		8 x UDR80-44	4 x 4	128	140 x 567 x 339
		Total: 16 racks	-	(=224)	-
<b>ULUF P610</b>	80	15 x UDR80-43	4 x 3	180	140 x 567 x 254
		5 x UDR80-44	4 x 4	80	140 x 567 x 339
		Total: 20 racks	-	(=260)	-
<b>ULUF P610 GG</b>	80	10 x UDR80-43	4 x 3	120	140 x 567 x 254
		10 x UDR80-44	4 x 4	160	140 x 567 x 339
		Total: 20 racks	-	(=280)	-
<b>ULUF P820</b>	80	21 x UDR80-43	4 x 3	252	140 x 567 x 254
		7 x UDR80-44	4 x 4	112	140 x 567 x 339
		Total: 28 racks	-	(=364)	-
<b>ULUF P820 GG</b>	80	14 x UDR80-43	4 x 3	168	140 x 567 x 254
		14 x UDR80-44	4 x 4	224	140 x 567 x 339
		Total: 28 racks	-	(=392)	-
<b>ULUF 15</b>	80	-	-	3	-
<b>ULUF 65</b>	80	2 x UDR80-33	3 x 3	18	140 x 434 x 254
<b>ULUF 60</b>	80	2 x UDR80-33	3 x 3	18	140 x 434 x 254
<b>ULUF 120</b>	80	4 x UDR80-33	3 x 3	36	140 x 434 x 254
<b>ULUF 125</b>	80	4 x UDR80-33	3 x 3	36	140 x 434 x 254
<b>ULUF 400</b>	80	15 x UDR80-43	4 x 3	180	140 x 567 x 254
<b>ULUF 450-2M</b>	80	15 x UDR80-43	4 x 3	180	140 x 567 x 254
<b>ULUF 500</b>	80	20 x UDR80-43	4 x 3	240	140 x 567 x 254
<b>ULUF 550-2M</b>	80	20 x UDR80-43	4 x 3	240	140 x 567 x 254
<b>ULUF 700</b>	80	25 x UDR80-43	4 x 3	300	140 x 567 x 254
<b>ULUF 750-2M</b>	80	25 x UDR80-43	4 x 3	200	140 x 567 x 254

Box size 90/96/100

Model	Box size (height mm)	Rack model & total quantity	No. of boxes per rack (depth x height)	Total no. of boxes	Exterior dimensions of rack (wxdxh mm)
<b>MUF 40</b>	90/96/100	2 x UDR100-22	2 x 2	8	140 x 294 x 212
<b>ULUF P10</b>	90/96/100	-	-	3	-
<b>ULUF P50</b>	90/96/100	2 x UDR100-32	3 x 2	12	140 x 420 x 212
<b>ULUF P90</b>	90/96/100	4 x UDR100-32	3 x 2	24	140 x 420 x 212
<b>ULUF P390</b>	90/96/100	9 x UDR100-42	4 x 2	72	140 x 558 x 212
		3 x UDR100-43	4 x 3	36	140 x 558 x 318
		Total: 12 racks	-	(=108)	-
<b>ULUF P390 GG</b>	90/96/100	6 x UDR100-42	4 x 2	48	140 x 558 x 212
		6 x UDR100-43	4 x 3	72	140 x 558 x 318
		Total: 12 racks	-	(=120)	-
<b>ULUF P500</b>	90/96/100	12 x UDR100-42	4 x 2	96	140 x 558 x 212
		4 x UDR100-43	4 x 3	48	140 x 558 x 318
		Total: 16 racks	-	(=144)	-
<b>ULUF P500 GG</b>	90/96/100	8 x UDR100-42	4 x 2	64	140 x 558 x 212
		8 x UDR100-43	4 x 3	96	140 x 558 x 318
		Total: 16 racks	-	(=160)	-
<b>ULUF P610</b>	90/96/100	15 x UDR100-42	4 x 2	120	140 x 558 x 212
		5 x UDR100-43	4 x 3	60	140 x 558 x 318
		Total: 20 racks	-	(=180)	-
<b>ULUF P610 GG</b>	90/96/100	10 x UDR100-42	4 x 2	80	140 x 558 x 212
		10 x UDR100-43	4 x 3	120	140 x 558 x 318
		Total: 20 racks	-	(=200)	-
<b>ULUF P820</b>	90/96/100	21 x UDR100-42	4 x 2	168	140 x 558 x 212
		7 x UDR100-43	4 x 3	84	140 x 558 x 318
		Total: 28 racks	-	(=252)	-
<b>ULUF P820 GG</b>	90/96/100	14 x UDR100-42	4 x 2	112	140 x 558 x 212
		14 x UDR100-43	4 x 3	168	140 x 558 x 318
		Total: 28 racks	-	(=280)	-
<b>ULUF 15</b>	90/96/100	-	-	3	-
<b>ULUF 65</b>	90/96/100	2 x UDR100-32	3 x 2	12	140 x 420 x 212
<b>ULUF 60</b>	90/96/100	2 x UDR100-32	3 x 2	12	140 x 420 x 212
<b>ULUF 120</b>	90/96/100	4 x UDR100-32	3 x 2	24	140 x 420 x 212
<b>ULUF 125</b>	90/96/100	4 x UDR100-32	3 x 2	24	140 x 420 x 212
<b>ULUF 400</b>	90/96/100	15 x UDR100-42	4 x 2	120	140 x 558 x 212
<b>ULUF 450-2M</b>	90/96/100	15 x UDR100-42	4 x 2	120	140 x 558 x 212
<b>ULUF 500</b>	90/96/100	20 x UDR100-42	4 x 2	160	140 x 558 x 212
<b>ULUF 550-2M</b>	90/96/100	20 x UDR100-42	4 x 2	160	140 x 558 x 212
<b>ULUF 700</b>	90/96/100	25 x UDR100-42	4 x 2	200	140 x 558 x 212
<b>ULUF 750-2M</b>	90/96/100	25 x UDR100-42	4 x 2	200	140 x 558 x 212